SISTEMA ECMOLIFE: DESCRITIVO TÉCNICO E MATERIAIS

REVISÃO 05 (12.05.2022)

SUMÁRIO

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REFERÊNCIAS

1. S. BOTTRELL et al. **A comparison study of haemolysis production in three contemporary centrifugal pump.** Perfusion 2014, Vol. 29 (5) 411-416.
2. CONRAD S, RYCUS P. **The registry of the extracorporeal life support organization.** In: Annich GL, Lynch WR, MacLaren G, Wilson JM, Barlett RH. (eds) ECMO. Extracorporeal cardiopulmonary support in critical care. 4th edition. 2012, pp 102-103.
3. RANJIT J. et al. **Outcomes of a multicenter trial of the Levitronix CentriMag ventricular assist system for short-term circulatory support.** The Journal of Thoracic and Cardiovascular Surgery, April 2011.
4. R. LORUSSO et al. **Effects of phosphorylcholine coating on extracorporeal circulation management and post-operative outcome: a double-blind randomized study.** Interactive Cardiovascular and Thoracic Surgery 8 (2009) 7-11.
5. DE SOMER F. FRANCOIS K, VAN OEVEREN W, POELAERT J, DE WOLF D. EBELS T, VAN NCOTEN G. **Phosphorylcholine coating of extracorporeal circuits provide natural protection against blood activation by the material surface.** Eur J. Cardio-Thorac 2000; 18: 602-606.
6. S. CALVI et al. **55 days ECMO support for post-cardiotomy, by-ventricular failure: a case report.** Abstract presented at 14th FECECT congress, 2011.
7. NAGI EL SHOUKI, M. STILLER. **Case report of 19 days extracorporeal life support with a new ECMO oxygenator.** Abstract presented at 14th FECECT congress, 2011.
8. M. PIERI et al. **A new phosphorylcholine-coated polymethylpentene oxygenator for extracorporeal membrane oxygenation: a preliminar experience.** Perfusion 2012; Volume: 28 issue: 2, page(s): 132-137.
9. F. PAPPALARDO et al. **Bioline® heparina-coated ECMO with bivalirudin anticoagulation in a patient with acute heparina-induced thrombocytopenia: the immune reaction appeared to continue unabated.** Perfusion 2009; 24; 135-137.
10. M. BELLATO et al. **A brief clinical case report of real-time monitoring of oxygenator performances and patient-machine interdependency during a prolonged VV ECMO support.** J Clin Monit Comput. 2017 Oct; 31 (5); 1027-1033.

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